

## Sóttgreinsisápa RBS 50



Rating: Not Rated Yet

[Fyrirspurn varðandi þessa vöru](#)

Framleiðandi: [RBS hreinsiefni](#)

### Vörulýsing

Fljótandi sápa með virkni á bakteríur og veirur fyrir áhöld sem lögð eru í bleyti, þvegin í þvottavél eða ultrasonísk böð. Hreinsar fljótt og vel blöð og vefjaleifar, öllu og flúu – jafnvel þó óhreinsindin hafi legið tengi á áhaldinu. Eftir þvott verða áhöldin glansandi og án nokkurra ráka. Mjög notadryggt, 3 ml af sápu í hvern lítra vatns. 5 Lítrar.

Sápan uppfyllir Evrópustaða EN 1040, EN 1276, EN 1275 og NF T 72-190 (spectrum 4) um hreinsun áhalds, glerlita o.s.frv. í iðnað og framleiðslu.

RBS 50 MD meets the requirements of the European Directive related to medical devices and is certified CE 1639.

- EN 1040 : basic bactericidal activity (Pseudomonas Aeruginosa, Staphylococcus Aureus)
  - EN 1276 : bactericidal activity in dirty conditions (3 g/L bovine albumine), (Pseudomonas Aeruginosa, Staphylococcus Aureus, Escherichia Coli, Enterococcus Hirae) [100 ml of RBS 50 MD in 1 l of water at 40°C – 60 minutes].
  - EN 1275 : basic fungicidal activity (Candida Albicans, Aspergillus Niger) [15 ml of RBS 50 MD in 1 l of water at 20°C - 60 min].
  - NFT 72-190 spectrum 4 : bactericidal activity [25 ml of RBS 50 MD in 1 l of water at 50°C - 15 min].
- A 2.5% (v/v) RBS 50 MD solution totally inhibits the infecting power of HIV virus.